# Digital Image Processing

From

Image processing

To

Deep-learning using CNN




IMAGE → PROCESSING → ANALSYSIS → RECOGNIZNG



AASTMT Smart Village

**Presented by: Mohamed Hazem** 

Supervised by: Engy Ehab

# Acknowledgement

I would like to thank our advisor Assoc. Prof. Dr Aliaa Youssif who really directed me to produce this project. I would like to give special thanks to our Professor for:

- Providing me with the needed technical information throughout the project.
- ❖ Not only the technical information but also for the soft skills; by giving me the chance to present every step of project in a formal way.
- ❖ Giving me the chance to search for the needed information in design.

Nonetheless, I would also thank Engineer Engy Ehab who gives me much knowledge in

- \* CNN (Convolutional neural network).
- ❖ Deep-learning as all
- Python
- Techniques

## Our Project Split to 3 Phases 2 of them Image Processing and Last one using CNN Deep-Learning.

No worry everything will be present clearly in Table of Contents

# Table of Contents

### Table of Contents

Chapter One	3
Analysis of MRI Images Using Image Processing Technique	
Example ONE	
Example Two	
Chapter Two	
Analysis of Breast Cancer Images Using Image Processing Technique	
Chapter Three & Last	
White Blood cells Classification using CNN Deen learning & TensorFlow & Keras & Rea	

## Chapter One


Analysis of MRI Images Using Image Processing Technique

#### **Example ONE**

So, here's our image and we going to use some process on it to get it out

Original Image





**Target Image** 



# Let's go ahead and see what need to happen Process & Pipeline:

#### **Step 1 Contrast:**

imgcon=cv2.addWeighted(imgabs,alpha,np.zeros(imgabs.shape, imgabs.dtype),0,beta)

**Note (Alpha control in Degree of Contrast)** 



#### Step 2 Blurring:

imggau2 = cv2.GaussianBlur(imgcon,(3,3),sigmaX=1,borderType=cv2.BORDER\_CONSTANT)



### **Step 3 Edge Detection using Canny:**

imgcanny= cv2.Canny(imggau2,300,800)



#### **Step 4 Repeat Step 3 with different values:**


After applying many of canny but with different values we got our Target Finally


As you can See →




#### **Example Two**

## We know What happens from Example one so we will move Quickly on this one 😊










Final Output



# **Chapter Two**

Analysis of Breast Cancer Images Using Image Processing Technique




At Target Image you can see the cancer clearly after Applying some Image Processing Techniques

### **Process & Pipeline:**

**Step 1 Negative:** 

imgneg = 255-img



### **Step 2 Dilation:**

imgdil = cv2.dilate(imgneg, kernel, iterations=1)




### Step 3 erosion:

imgero = cv2.erode(imgdil, kernel, iterations=1)



### **Step 4 Threshold:**

imgthreshold = cv2.threshold(imgero,50,90,cv2.THRESH\_BINARY)



## **Chapter Three & Last**

White Blood cells Classification using CNN Deep learning & TensorFlow & Keras & Real-life Dataset

First, Let's Define our Dataset

It's 9956 Real images for white Blood Cells (EOSINOPHIL, LYMPHOCYTE, MONOCYTE, NEUTROPHIL)


We now will pick some of it to train and other to test and also a few as validation


### **Train**

Train\_Data\_Path = Path("C:/Users/moham/Downloads/dataset2-master/dataset2-master/images/TRAIN")

```
print("EOSINOPHIL: ", Train_JPG_Labels.count("EOSINOPHIL"))
print("LYMPHOCYTE: ", Train_JPG_Labels.count("LYMPHOCYTE"))
print("MONOCYTE: ", Train_JPG_Labels.count("MONOCYTE"))
print("NEUTROPHIL: ", Train_JPG_Labels.count("NEUTROPHIL"))

EOSINOPHIL: 2497
LYMPHOCYTE: 2483
MONOCYTE: 2478
NEUTROPHIL: 2499
```







#### **Test**

Test\_Data\_Path = Path("C:/Users/moham/Downloads/dataset2-master/dataset2-master/images/TEST")

```
print("EOSINOPHIL: ", Test_JPG_Labels.count("EOSINOPHIL"))
print("LYMPHOCYTE: ", Test_JPG_Labels.count("LYMPHOCYTE"))
print("MONOCYTE: ", Test_JPG_Labels.count("MONOCYTE"))
print("NEUTROPHIL: ", Test_JPG_Labels.count("NEUTROPHIL"))
```

EOSINOPHIL: 623 LYMPHOCYTE: 620 MONOCYTE: 620 NEUTROPHIL: 624







## Validation

Validation\_Data\_Path = Path("C:/Users/moham/Downloads/dataset2-master/dataset2-master/images/TEST\_SIMPLE")


```
print("EOSINOPHIL: ", Validation_JPG_Labels.count("EOSINOPHIL"))
print("LYMPHOCYTE: ", Validation_JPG_Labels.count("LYMPHOCYTE"))
print("MONOCYTE: ", Validation_JPG_Labels.count("MONOCYTE"))
print("NEUTROPHIL: ", Validation_JPG_Labels.count("NEUTROPHIL"))

EOSINOPHIL: 13
LYMPHOCYTE: 6
MONOCYTE: 4
NEUTROPHIL: 48
```







# Let's take a look on Our images that we will train



# And Our images that we will test



# **Here also Our Validation images**



### Now we will class Our Data as numbers

```
TRAIN:
{'EOSINOPHIL': 0, 'LYMPHOCYTE': 1, 'MONOCYTE': 2, 'NEUTROPHIL': 3}
[3, 2, 0, 0, 3]
(220, 220, 3)

VALIDATION:
{'EOSINOPHIL': 0, 'LYMPHOCYTE': 1, 'MONOCYTE': 2, 'NEUTROPHIL': 3}
[3, 3, 3, 3, 1]
(220, 200, 3)

TEST:
{'EOSINOPHIL': 0, 'LYMPHOCYTE': 1, 'MONOCYTE': 2, 'NEUTROPHIL': 3}
[3, 1, 1, 1, 0]
(220, 220, 3)
```

So, we can train it and get result as Predictions of (0,1,2,3)

Though...

Things will Happens at next pages

## Let's the Deep-learning Start

### **Activation of CNN**

### \*CNN STRUCTURE WITH SeparableConv2D\*

```
Model = Sequential()
Model.add(SeparableConv2D(32,3,
                          activation="relu",
                 input_shape=(220,220,3)))
Model.add(BatchNormalization())
Model.add(MaxPooling2D((2)))
Model.add(SeparableConv2D(64,3,
                activation="relu"))
Model.add(SeparableConv2D(128,(3,3),
                activation="relu"))
Model.add(Dropout(0.5))
Model.add(MaxPooling2D((2)))
Model.add(SeparableConv2D(64,3,
                activation="relu"))
Model.add(SeparableConv2D(128,3,
                activation="relu"))
Model.add(Dropout(0.5))
Model.add(GlobalAveragePooling2D())
Model.add(Flatten())
Model.add(Dense(256,
               activation="relu"))
Model.add(Dropout(0.5))
Model.add(Dense(4,
                activation="softmax"))
```

### Then we add our TensorFlow Keras

```
Call_Back = tf.keras.callbacks.EarlyStopping(monitor="val_accuracy",patience=5,mode="max")
```

### And here just to define loss & accuracy print in our cmd

```
{\tt Model.compile} (optimizer = "rmsprop", loss = "categorical\_crossentropy", metrics = ["accuracy"])
```

## So, Let's train it up

# Normal we use between 10 epochs to 28 epochs But to get best result & accuracy we will go for 50 epochs

## Small look on training...

```
Epoch 1/50
312/312 [===
   0.1714
Epoch 3/50
312/312 [============] - 645s 2s/step - loss: 0.6806 - accuracy: 0.7265 - val_loss: 0.7950 - val_accuracy: 0.
6000
Epoch 4/50
7429
Epoch 5/50
8571
Epoch 6/50
312/312 [===
   8286
```

# So, after 8 hours of training, we finally back......... Let's check accuracy of our model now


# And here's model summary if anyone interest

<pre>print(Model.summary())</pre>		
Model: "sequential"		
Layer (type)	Output Shape	Param #
separable_conv2d (Separable Conv2D)		
batch_normalization (BatchN ormalization)	(None, 218, 218, 32)	128
<pre>max_pooling2d (MaxPooling2D )</pre>	(None, 109, 109, 32)	0
separable_conv2d_1 (Separab leConv2D)	(None, 107, 107, 64)	2400
separable_conv2d_2 (Separab leConv2D)	(None, 105, 105, 128)	8896
dropout (Dropout)	(None, 105, 105, 128)	0
<pre>max_pooling2d_1 (MaxPooling 2D)</pre>	(None, 52, 52, 128)	0
separable_conv2d_3 (Separab leConv2D)	(None, 50, 50, 64)	9408
separable_conv2d_4 (Separab leConv2D)	(None, 48, 48, 128)	8896
dropout_1 (Dropout)	(None, 48, 48, 128)	0
global_average_pooling2d (G lobalAveragePooling2D)	(None, 128)	0
flatten (Flatten)	(None, 128)	0
dense (Dense)	(None, 256)	33024
dropout_2 (Dropout)	(None, 256)	0
dense_1 (Dense)	(None, 4)	1028
Total params: 63,935 Trainable params: 63,871 Non-trainable params: 64		


# So, let's start now with some of CNN Benefits

## Note (we got best result at 23 epochs so it didn't continue)


# **Accuracy Plot**




## **Loss Plot**



## Val Acc vs Val Loss



# And here's everything




So, let's go ahead and have some Image Processing

## Let's Pick a random Image and Play with it

```
Any_IMG = Main_Train_Data["JPG"][6]
IMG = image.load_img(Any_IMG,target_size=(220,220))
Array_IMG = image.img_to_array(IMG)
Array_IMG = np.expand_dims(Array_IMG,axis=0)
Array_IMG /= 255


plt.imshow(Array_IMG[0])
plt.title("ANY TARGET IMAGE")
plt.show()
```





## **Apply some layers?**


# We can apply 0-15 different layer


### Here's some





























So, we Almost at the End 😕

If you remember we Classed our White Blood Cells as 0,1,2,3

# So, here's the Predication & Classification now

(EOSINOPHIL 0, LYMPHOCYTE 1, MONOCYTE 2, NEUTROPHIL 3)

